Introducing Numera Libris #1 In Health & Safety Wearables





SAN ANTONIO AUSTIN 210-807-8488 / 512-402-8999

CORPUS CHRISTI 361-808-4828



About the Device

- Portable On The Go Safety-Completely mobile on the AT&T Cellular Network with GPS Tracking
- Advanced Fall Detection-Designed to alert a central station even if you are not able to.
- Never Alone-Instant two-way hands-free voice communication allows the individual to speak to the monitoring center directly from the portable device.
- Use Daily-Water Resistant up to 3 feet for 30 minutes (IP67 rating); Shatter Resistant from 6 feet (IK06 rating)
- 36-Hour Battery Life

numera Libris™






Health & Safety Management Users


- Seniors & Assisted Living
- Construction Workers
- Realtors
- Legal & Medical Professionals
- Service & Field Technicians
- Teachers & Administrators
- Law Enforcement
- Everyone



Features & Functions



- Cellular Indicator shows cellular signal strength
- 2. Battery Indicator, which surrounds the Call Button, shows battery status
- 3. Call Button, when pushed, initiates a two-way call to your response team
- 4. Speaker allows you to hear your response team



- 5. Microphone allows the response team to hear you
- 6.
- 7. Clip Lock enables you to affix a Lanyard/Clip or Belt Clip
- 8.
- 9. Charging Pins connect Libris to the Charger



Large Predictive Fall Market

- 33% of adults 65+ fall each year
- 43% of falls occur outside the home or office
- After 1st fall, 3X more likely to fall again
- 40 million in '10; 55 million in '20
- 89% have cell phones; 9-13% are smartphone users
- 28 million households; 11 million live alone; 8 million women
- 96% want to live on their own, 1/3 more women than men

Sources: CDC



Smart Fall Detection Tested and Approved

Models built from 10,000 real human falls, not computer

simulations.

Has 4 sensors and

87% of all currently active

devices have Fall detection

enabled.

Cloud based analysis


provides limitless computing.

Our model is easily updated.



Ensuring a "fall" is not a "drop"

Data from monitored activity





Numera Libris Kit







Libris Charger











Power Cord

User Guide



Quick set-up simple steps



1 Press and hold the Call Button with the device in the Charger.




2 The response team will answer, and will walk the user through customizing Libris and completing set up.



3 Charge Libris



Processing an "event"



Libris detects suspicious activity

Sends data via cellular network

Analyzed by cloud-based classifier

If FALL, monitoring center receives alert... two-way talk to VERIFY

Monitoring center dispatches assistance, if needed

Contact List is notified



35-40 SECONDS







Equipment and Monthly Monitoring Rate

Equipment Cost

1 Numera Kit: \$200
Includes: Numera Libris, Charging

Station

Belt Clip and Lanyard

Monitoring Rate

- \$48.00 per month
 - Includes Cell and Data plan
 - GPS Tracking & Fall Detection
 - Central Station Monitoring
 - Access to the Numera EverThere Website Portal E-mail notification alerts






Introducing Numera EverThere Web Portal

Numera EverThere features include:

- Location of Device (if misplaced will locate by chime)
- Breadcrumbs of last 25 events
- Signal and Battery Strength
- Unified Log in with other Web services
- Volume adjustment
- Browser based





Summary

- Health and Safety Management
- Outperforms traditional PERS offerings: not just a panic button
- The most advanced fall detection available
- Central Station Monitored
- Built for growth
- Visit www.ender.services/numera

numera Libris™



